

DEVELOPMENT OF HYDROXYLATED POLYCYCLIC AROMATIC HYDROCERBON(OH-PAHs) ANALYTICAL METHOD IN BIOLOGICAL SAMPLE WITH UPLC/Q-TOF MS

Sato N^{1*}, Yonekubo J¹, Ezaki T¹, Suzuki M², Matsumura C³, Haga Y³, Nakano T^{3,4}

¹ Nihon Waters K.K., 1-3-12 Kitashinagawa, Shinagawa-ku Tokyo-to 140-0001 Japan;

² Kobe University, Graduate School of Agricultural Science. 1-1 Rokkodaicho, Nada, Kobe, Hyogo 657-8501, Japan;

³ Hyogo Prefectural Institute of Environmental Sciences, 3-1-27, Yukihiro-cho, Suma-ku, Kobe-city Hyogo-ken, 654-0037 Japan;

⁴ School of Engineering, Osaka University, Japan

Introduction

Polycyclic aromatic hydrocarbons (PAHs) are known as environmental contaminants which occur in oil, coal and tar deposits, and are produced as byproducts of fuel burning. PAHs are also found in cooked foods. Some of PAHs may cause abnormal effect such as carcinogenicity, mutagenicity and teratogenicity. Some studies determined the residue levels of PAHs and metabolized PAHs (hydroxylated PAHs : OH-PAHs) in biological sample. Hydroxyl group of OH-PAHs has high polarity and boiling point, therefore OH-PAHs was made derivatization such as acetylation, trimethylsilylation and analyzed with GC-MS. Complicated preparations are needed for analyzing with GC-HRMS, beside that assining metabolized PAHs and derivatized OH-PAHs is difficult.

For solving analytical issue of measurement of metabolized environmental contaminants, our previous study determined elution order for 51 congeners of OH-PCB without derivatization and developed an analytical method for quantity to separate mixture of 6 major OH-PCBs in human blood, then applied analytical method to biological sample of human urine and blood with UPLC/QTOF MS. As compound that might be hydroxypyrene was detected in urine sample, other OH-PAHs might be exist in biological sample. The present study aimed at developing analytical method of OH-PAHs that be able to separate each isomer, applying method to biological sample and measuring OH-PAHs for comprehensive evaluation of accumulation of environmental contaminants.

Materials and methods

OH-PAHs Standard Solution:

Ten standard solution of OH-PAHs (see Fig.1) were purchased from Cambridge Isotope Laboratories. Inc. and solution was diluted in acetone.




Fig.1 Ten compounds of OH-PAH

LC-TOF analysis:

Identification and quantification were performed using ultra performance liquid chromatography (UPLC: Waters Acquity UPLC system) and a high-resolution q-tof mass spectrometer (Xevo G2-S QToF MS) with a resolving power of more than 32500 (Table 1)

Provide a succinct description of the materials and methods used in your research².

Table 1 Analytical conditions of UPLC/Q-ToF MS

UPLC			MS	
Instrument:	ACQUITY UPLC		Instrument:	Xevo G2-S Q-TOF
Column:	BEH C18 2.1ID X 150 mm, 1.7um		Ionization mode:	ESI negative
Flow rate:	0.5 mL/min.		Capillary:	1.5 kV
Column heater:	60 degree celsius		Sampling Cone:	40 V
Mobile Phase A:	5mM CH ₃ COONH ₄ aq.		Source Temp:	120 degree celsius
Mobile Phase B:	THF/CH ₃ CN (v/v: 1/4)		Desolvation Temp:	600 degree celsius
Gradient :			Cone Gas Flow:	20 L/hr.
Time(min)	%A	%B	Desolvation Gas Flow:	800 L/hr.
0.0	75	25	Resolving Power:	<32,500
17.0	25	75	Selected Ion:	
18.0	1	99	Hydroxynaphthalene	<i>m/z</i> = 143.0497
18.5	75	25	Hydroxyfluorene	<i>m/z</i> = 181.0653
Total Run Time:	20 min		Hydroxyphenanthrene	<i>m/z</i> = 193.0653
			Hydroxypyrene	<i>m/z</i> = 217.0653

Results and discussion

This method applied for human biological samples. Mass chromatogram and spectra of compound that might be hydroxypyrene or hydroxyfluoranthene in human urine sample are shown in Fig.1 and Fig. 2.




Fig.1 Mass chromatogram of $m/z = 217.0653$ in urine sample




Fig.2 MS spectra of retention time 7.92 min in urine sample

Result of elemental composition of $m/z=217$ in urine sample is shown in Fig. 3. Composition was decided by accurate mass and isotope pattern. There were two candidates gratifying number of element, carbon, hydrogen, nitrogen and oxygen. This compound was decided C16H10O same as composition of hydroxypyrene and hydroxyfluoranthene by accurate mass and isotope pattern.




Fig.3 Result of elemental composition of spectra in urine sample at retention time of 7.92 min

References

- Hirai T., Kinoshita H., Okamura H., Yano Y., Nakano T., 2012 Development of Simultaneous Determination Method of Hydroxylated Polycyclic Aromatic Hydrocarbons in Urine by LC/MS/MS and Its Application to Assessment of Polycyclic Aromatic Hydrocarbons Exposure. BUNSEKI KAGAKU Vol. 61, No. 11, pp. 925-930.
- Boffetta P., Jourenkova N., Gustavsson P., 1997 Cancer risk from occupational and environmental exposure to polycyclic aromatic hydrocarbons. Cancer Causes & Control 1997, Volume 8, Issue 3, 444-472
- Grete J., Jonny B., David W., Freck A., 2003 The application of HPLC-F and GC-MS to the analysis of selected hydroxy polycyclic hydrocarbons in two certified fish bile reference materials. J. Environ. Monit., 2003,5, 513-520
- Sasaki H., Yonekubo J., Kanai M., Toriba A., Kizu R., Hayakawa K., 2002 Simultaneous detection of monohydroxybenzo(a)pyrene positional isomers by reversed-phase liquid chromatography coupled to electrospray ionization mass spectrometry. BIOMEDICAL CHROMATOGRAPHY, 16, 432-436
- Sato N., Yonekubo J., Ezaki T., Suzuki M., Matsumura C., Haga Y., Nakano T., 2012 Measurement of accumulation of hydroxylated polychlorinated biphenyl (OH-PCBs) in human urine and blood. Dioxin 2012
- Suzuki M., Okuno T., Matsumura C., Sato N., Yonekubo J., Ezaki T., Inoue Y., Imaishi H., Nakano T., 2011 Development of hydroxylated polychlorinated biphenyl(OH-PCB) analytical method in human urine with UPLC/Q-ToF MS. Dioxin 2011
- Kuroki H., Haraguchi K., Saito H., Masuda Y., Wehler EK., Bergman A., 1993 Accumulation of hydroxylated PCB metabolites in blood. Fukuoka Igaku Zasshi. May;84(5), 248-56